



USB-CAN-LAN-RS232

CONVERTER module

User's Manual

Rev. 1805

CONTENTS

| | | |
|------------------|---|-----------|
| CHAPTER 1 | GENERAL INFORMATION | 1 |
| 1.1. | SHORT DESCRIPTION | 1 |
| 1.2. | COMMUNICATION PROTOCOLS | 2 |
| 1.2.1. | <i>Ethernet (LAN, Wireless LAN (WLAN))</i> | 2 |
| 1.2.2. | <i>USB (USB DEVICE connector)</i> | 2 |
| 1.2.3. | <i>RS232 (RS232 connector)</i> | 2 |
| CHAPTER 2 | INSTALLATION | 3 |
| 2.1. | WINDOWS 7/8 | 3 |
| 2.2. | WINDOWS XP | 4 |
| CHAPTER 3 | REST API OVER LAN OR WLAN PROTOCOL | 5 |
| 3.1. | COMMUNICATION CONTROL | 5 |
| 3.2. | ADDRESSING CONVENTIONS | 5 |
| 3.3. | MODULES AND REGISTERS | 6 |
| 3.4. | PROTOCOL DESCRIPTION | 7 |
| 3.4.1. | <i>Commands and Queries</i> | 7 |
| CHAPTER 4 | ASCII SERIAL PROTOCOL | 9 |
| 4.1. | PHYSICAL CONNECTION | 9 |
| 4.2. | REGISTERS | 9 |
| 4.3. | PROTOCOL DESCRIPTION | 10 |
| 4.3.1. | <i>Commands and Queries</i> | 11 |
| CHAPTER 5 | SETUP OF LASER CONTROL OVER LAN OR WLAN | 13 |
| 5.1. | LAN CONFIGURATION | 13 |
| 5.2. | IMPORTANT NOTES | 14 |
| CHAPTER 6 | COMMAND ERRORS | 15 |
| CHAPTER 7 | COMMUNICATION TROUBLESHOOTING GUIDE | 17 |
| 7.1. | PHYSICAL CONNECTION | 17 |
| 7.2. | TESTING AND TROUBLESHOOTING PROCESS FLOW | 17 |
| 7.3. | COMMUNICATION TESTING UTILITY | 17 |
| 7.4. | STEP-BY-STEP INSTRUCTIONS | 17 |
| 7.5. | TESTING LAN WITHOUT AN USB/RS232 CONNECTION | 19 |
| CHAPTER 8 | TABLES | 21 |
| 8.1. | REST API OVER LAN OR WLAN REGISTER LIST EXAMPLE | 21 |
| 8.2. | ASCII SERIAL PROTOCOL REGISTER LIST EXAMPLE | 22 |
| 8.3. | LIST OF COMMANDS FOR REST OVER API PROTOCOL | 23 |
| 8.4. | LIST OF COMMANDS FOR ASCII SERIAL PROTOCOL | 29 |

LIST OF FIGURES

| | |
|---|----|
| FIGURE 1 LASER/BLOCK HIERARCHIES FOR MULTI-SERVER (LEFT) AND STAND-ALONE (RIGH) SYSTEMS | 5 |
| FIGURE 2 RS232 CONNECTOR PINOUT | 9 |
| FIGURE 3 MAIN SETUP PAGE | 13 |
| FIGURE 4 <i>WORKING DIRECT RS232 INTERFACES FOUND LIST</i> | 18 |
| FIGURE 5 SUCCESSFUL CONNECTION RESPONSE | 18 |
| FIGURE 6 RS232 CONNECTOR PINOUT ON THE LASER SIDE | 19 |
| FIGURE 7 LAN SCAN RESULTS | 20 |

LIST OF TABLES

| | |
|---|----|
| TABLE 1 COLUMNS OF THE REGISTER LIST FOR REST API | 6 |
| TABLE 2 REGISTER DATA TYPES FOR <i>TYPE</i> COLUMN | 7 |
| TABLE 3 COMMANDS AND QUERIES FOR REST API | 8 |
| TABLE 4 COLUMNS OF THE REGISTER LIST FOR ASCII | 10 |
| TABLE 5 COMMANDS AND QUERIES FOR ASCII | 11 |
| TABLE 6 REST API OVER LAN OR WLAN REGISTER LIST EXAMPLE | 21 |
| TABLE 7 ASCII SERIAL PROTOCOL REGISTER LIST EXAMPLE | 22 |
| TABLE 8 LIST OF COMMANDS FOR REST OVER API PROTOCOL | 23 |
| TABLE 9 LIST OF COMMANDS FOR ASCII SERIAL PROTOCOL | 29 |

This document describes remote control of EKSPLA products using the USB-CAN-LAN-RS232 CONVERTER Module.

The Module is installed in EKSPLA products either as an optional add-on or as a standard feature. If the Module is provided as an optional add-on, its interfaces will function in parallel with the native USB/RS232 interfaces of the laser. Native USB/RS232 interfaces use different drivers and protocols and should be not confused with Module interfaces.

This document covers interfaces associated with the Converter module only.

1.1. Short Description

The Module acts as a communication bridge allowing control of various Ekspla products over a number of hardware interfaces:

- RS232
- USB
- Local Area Network (LAN)
- Wireless LAN (WLAN)

In addition to communication functions, the Module can log operational parameters to an SD card or USB stick.

The Module is a single printed circuit board (PCB), mounted either into the laser power supply unit (PSU), or into a separate enclosure:

- Internally mounted modules are installed into many types of laser PSUs, e.g. PS8000, PS8600, PS6000. The Module is distinguished from other interfaces by connectors named USB HOST, USB DEVICE, LAN and RS232.
- The standalone Converter module type is named 'USB-CAN-LAN-RS232 CONVERTER'.

A major advantage of using the Converter is that it allows controlling and diagnosing the laser from virtually any device: a microprocessor, PC, a smartphone, running any operating system (Windows, iOS, Linux etc.) No software driver installation is necessary in most cases, except USB connection to Windows versions XP...8 (some exceptions apply).

No special application is needed to control the laser. Because the communication is text-based, standard applications, such as 'terminal' for USB/RS232 or a web browser for LAN, are sufficient.

The Converter module can bridge the following hardware interfaces:

- USB virtual serial port – EKSPLA CAN bus
- RS232 serial port – EKSPLA CAN bus

- 10BASE-T/100BASE-TX based LAN – EKSPLA CAN bus

The EKSPLA CAN bus is an internal bus that interconnects all laser drivers and sensors. The laser system can be controlled and diagnosed by sending messages via this bus. Although direct bridges work only in direction from/to the EKSPLA CAN bus, the messages on the bus are visible from all interfaces. This is useful for communication debugging, because it allows bus sniffing via a different interface.

1.2. Communication Protocols

1.2.1. Ethernet (LAN, Wireless LAN (WLAN))

- REST API running over HTTP protocol (HTTP GET command).

HTML is used for messages; therefore any web browser may be used for communication testing.

- CAN messages tunneling over TCP/IP

Closed protocol. *REMOTECONTROL.dll* is used on the PC side to communicate with the laser. The *CAN Browser* application uses this mode of communication.

1.2.2. USB (USB DEVICE connector)

- Protocol based on ASCII readable strings

Plain text messages are used for messaging; therefore, any terminal application is suitable for communication testing.

- CAN messages tunneling over a USB virtual serial port

Closed protocol. *REMOTECONTROL.dll* is used on the PC side to communicate with the laser.

1.2.3. RS232 (RS232 connector)

- Protocol based on ASCII readable strings

Plain text messages are used for messaging; therefore, any terminal application is suitable for communication testing.

- CAN messages tunneling over a USB virtual serial port

Closed protocol. *REMOTECONTROL.dll* is used on the PC side to communicate with the laser.

USB and RS232 share the same protocols.

Note:

Due to limited support for MS Windows 7 and earlier versions, driver installation and proper operation is guaranteed to work only with MS Windows 10.

No special installation procedures are needed for:

- LAN, Wireless LAN and RS232 – on all systems.
- USB – on MS Windows 10, Linux.

USB on Windows XP, Windows Vista, Windows 7 and Windows 8 requires installation of a USB-serial driver. This happens automatically for some computers connected to the Internet, as soon as the USB cable from the Module is connected to the PC. In case this automatic installation fails, a manual procedure will have to be performed.

First of all, download the required file; click the following link, select 'Save as', save the file to disk:

<https://www.kernel.org/doc/Documentation/usb/linux-cdc-acm.inf>

It is also possible to enter the link into a web browser, copy the contents of the displayed page to a blank *Notepad* document and save it as *linux-cdc-acm.inf* (make sure to save it in the *.inf* format).

2.1. Windows 7/8

1. Open the *Start menu* and:
 - a. Windows 7: open *Control Panel* → *System and Security* → *System* → *Device Manager*.
 - b. Windows 8: point to *Settings*, select *Control Panel*, select *System*. Then, on the left, click *Device Manager*.
2. In *Ports (COM & LPT)*, you should see an open port *Gadget Serial V2.4*.

Please note the number of this port – it is relevant when connecting the *CAN Browser* application to your product.
3. Right-click *Gadget Serial V2.4*, choose the *Update Driver Software* option.
4. Select *Browse my computer for driver software*.
5. Navigate to the location of the downloaded file by selecting the *Let me pick from a list of device drivers on my computer* option, then *Show All Devices*, then *Have disk...*
6. Windows will display a message that the driver did not pass Windows Logo testing; ignore it, select *Continue anyway* and finish the driver installation.

2.2. Windows XP

1. After plugging in the module, the *Found New Hardware Wizard* starts up.
2. Select *Install from a list or specific location (Advanced)*.
3. Select *Include this location in the search* and enter the path or browse to the folder containing the file.
4. Windows will display a message that the driver did not pass Windows Logo testing; select *Continue anyway* and finish the driver installation.
5. Go to *Start Menu → Control panel → System → Hardware → Device Manager*.
6. In the *Device Manager*, expand *Ports (COM & LPT)* – *Gadget Serial* should be listed as the driver for one of the COM ports.

Once the drivers have been successfully installed, the USB to PC connection will appear as a COM port.

Use the LAN connector on the front panel of the Module to physically connect to the Module. The WiFi module is optional.

See Chapter 5 for setup instructions.

3.1. Communication Control

The communication is controlled by a LAN/WLAN/USB/RS232/CAN communication server. It links different buses by relaying messages, keeps the system configuration file(s), and stores messages in queues to enable data transfer between buses with different speeds and latency times. More complex systems may have multiple servers attached to several CAN bus segments.

Figure 1 illustrates how individual entities are addressed within a multi-server (left) and a stand-alone (right) system.




Figure 1 Laser/block hierarchies for multi-server (left) and stand-alone (right) systems

3.2. Addressing Conventions

Servers are addressed either by IP or by host name. REST requests are served by a special low-latency HTTP server on the 8080 port. Internally, the application software uses *Device ID* for block/laser identification.

Modules are addressed by name and ID in the form */(Module name)/(ID)*.

Registers are addressed by name.

3.3. Modules and Registers

The internal block/laser hierarchy is stored in a special configuration file REMOTECONTROL.CSV. The CAN Browser service utility reads the hierarchy on the fly, the rest of the software uses the stored REMOTECONTROL.CSV file.

Table 1 describes the columns of the register list.

Table 2 lists the data types for the *Type* column of the register list.

Table 6 (see Chapter 8) is an example of a laser register list. The actual register list of your product can be found in the software CD folder `\CAN network\Remote Control test utility\REMOTECONTROL.csv` or can be downloaded from the server over the LAN Connection at <http://your host or ip/REMOTECONTROL.CSV>.

Table 1 Columns of the register list for REST API

| Column | Description |
|-----------------------|--|
| Name/ID | Each module has an ID and a name. ID is a number in the range 0...63. (decimal), but in the table, it is given in HEX form, e.g. \$10 in HEX is 16 in DEC. In commands, the decimal value of ID is used, and it is tied to the name. If the product is comprised of a single module only, the command will always begin with the same symbols: /LDD1A/18. |
| Menu | Not used in remote control. |
| Type | Type of data in the register, see Table 2. Protocol operates in numbers converted to floating point and text representation; therefore the type may be ignored. |
| User rights | In remote control, value <i>ArUrSr</i> marks a read-only parameter. All other values are for writable parameters. |
| Non-volatile | The <i>NV</i> value indicates that a parameter may be written to non-volatile memory. |
| Min value / Max value | Parameter bounds. If exceeded, the communication controller will respond with an error message. |
| Short name | Not used. |
| Print format | Integer. <i>%u[x]</i> , where x is the dimension of parameter. Decimal. Because the register value is an integer, for display purposes it must be divided by a factor, e.g. <i>%.1f[x]</i> – by 10, <i>%.2f[x]</i> – by 100, etc. Float. <i>%f[x]</i> Set type [XXX, YYY, ...], with register type <i>u8</i> . The register value points to the actual (n-1) th element of the set. In ASCII remote control, actual set elements are sent, not the register value, e.g. if the print format is [OFF, ON], the reply to an inquiry will be 'OFF' or 'ON', not 0 or 1. |
| Name | Full register name, used to address the parameter. Copy this value to <i>Register name</i> argument when making function calls. |
| Value | Value of register. |

Table 2 Register data types for *Type* column

| Type | Description |
|-------------|--|
| u8 | byte, 0..255, unsigned 8-bit |
| s8 | shortint, -128..127, signed 8-bit |
| u16 | word, 0..65535, unsigned 16-bit |
| s16 | smallint, -32768..32767, signed 16-bit |
| u32 | longword, 0..4294967295, unsigned 32-bit |
| s32 | long integer, -2147483648..2147483647, signed 32-bit |
| float | single, 1.5×10^{-45} .. 3.4×10^{38} , Significant digits 7-8 |
| string8 | array of bytes |

3.4. Protocol Description

Commands and queries are sent by issuing HTTP GET requests, which are not authenticated.

The command line is an URL string.

The response body is an html document.

URL starts with the host name (IP address) and port, e.g. <http://host:8080> or <http://ip:8080>. The host name is programmable. By default, it is the product serial number. Please note that not all networks resolve the host name properly. See Chapter 5 for details about host names and IP addresses.

All of the following examples are shown without the host address part, e.g. instead of [http://192.168.1.192:8080/list\(\)](http://192.168.1.192:8080/list()), you will see just /list().

Command/query groups:

- Server commands/queries - /ID(), /List(), /Logget(), /Logstart()

These are not relayed to the CAN bus and work even with no modules live on the CAN bus.

- CAN module commands/queries - /(Module name)/(ID)/(Register name)

Modules are commanded and queried by performing register read and write operations.

3.4.1. Commands and Queries

A short list of commands and queries is presented in Table 3.

See Table 8 in Chapter 8 for an expanded list of commands with their sample responses. See Chapter 6 Command Errors for a list of command errors.

Table 3 Commands and queries for REST API

| Name | Syntax |
|--|---|
| Communication test | / |
| Device ID | /id() |
| Register list query | /list() |
| Read command | /(Module name)/(ID)/(Register name) |
| Write command | /(Module name)/(ID)/(Register name)/(Value) |
| Write to NVRAM command | /(Module name)/(ID)/(Register name)/(Value)/NV |
| Start/stop logging register updates | /(Module name)/(ID)/(Register name)/logstart()/(buffer size in bytes) Here: (buffer size in bytes) > 0 will turn on logging register changes (buffer size in bytes) = 0 will turn off logging register changes |
| Receive values with timestamp from the log | /(Module name)/(ID)/(Register name)/logget()/(number of records) |

This chapter describes remote control of Ekspla products via RS232 or USB DEVICE ports.

The commands shown here are just examples. For relevant commands, see the REMOTECONTROL.CSV configuration file of your product.

4.1. Physical Connection

To connect, use the RS232 female type connector or USB DEVICE connector on the front panel of your product.

The parameters of RS232 port are as follows:

- Speed: 19200
- Data bits: 8
- Stop bits: 1
- Parity: none
- Flow control: none

Figure 2 shows the RS232 connector pinout between the PC and the controller. The cable is terminated by a male-female D-Sub 9 pin connector.




Figure 2 RS232 connector pinout

4.2. Registers

The internal block/laser hierarchy is stored in a special REMOTECONTROL.CSV configuration file. The CAN Browser service utility reads the hierarchy on the fly, the rest of the software uses the stored REMOTECONTROL.CSV file.

Table 3 is an example of a laser register list. The actual register list of your product can be found in the software CD folder `\CAN network\Remote Control test utility\REMOTECONTROL.csv` or can be downloaded from the server over the LAN Connection at <http://your host or ip/REMOTECONTROL.CSV>.

Table 4 describes the register columns.

Table 7 (see Chapter 8) is an example of a laser register list

Table 4 Columns of the register list for ASCII

| Column | Description |
|-------------------------|--|
| Module name / Module ID | Each module has an ID and a name. ID is a number in the range 0...63 (decimal). Module name and ID are tied together in commands, e.g. /SY3PL50M/32. |
| Type | Type of data in the register, see Table 2. Protocol operates in numbers converted to floating point and text representation; therefore the type may be ignored. |
| User rights | In remote control, value <i>ArUrSr</i> marks a read-only parameter. All other values are for writable parameters. |
| Non-volatile | The <i>NV</i> value indicates that a parameter may be written to non-volatile memory. |
| Min value / Max value | Parameter bounds. If exceeded, the communication controller will respond with an error message. |
| Print format | Integer. <i>%u[x]</i> , where x is the dimension of parameter. Decimal. Because the register value is an integer, for display purposes it must be divided by a factor, e.g. <i>%.1f[x]</i> – by 10, <i>%.2f[x]</i> – by 100, etc. Float. <i>%f[x]</i> Set type [XXX, YYY, ...], with register type <i>u8</i> . The register value points to the actual (n-1) th element of the set. In ASCII remote control, actual set elements are sent, not the register value, e.g. if the print format is [OFF, ON], the reply to an inquiry will be 'OFF' or 'ON', not 0 or 1. |
| Register name | Full register name, used to address the parameter. Copy this value to <i>Register name</i> argument when making function calls. |
| Captured value | Value of register. |
| Comments | Comments. |

4.3. Protocol Description

The command is an ASCII string.

Symbols used in communication:

- #13 – carriage return [CR] ('Enter' key)
- #10 – line feed
- #03 – end of text

The command is entered by a #13 [CR] symbol.

#13#10 signals the end of line.

#03 ends the message.

Three apostrophes in a row (' ' ') precede a command error message.

4.3.1. Commands and Queries

A short list of commands and queries is presented in Table 5. See Table 9 in Chapter 8 for an expanded list of commands with their sample responses. See Chapter 6 Command Errors for a list of command errors.

Table 5 Commands and queries for ASCII

| Name | Syntax |
|--|---|
| Communication test | / |
| Device ID | /id() |
| Register list query | /list() |
| Read command | /(Module name)/(ID)/(Register name) |
| Write command | /(Module name)/(ID)/(Register name)/(Value) |
| Write to NVRAM command | /(Module name)/(ID)/(Register name)/(Value)/NV |
| Start/stop logging register updates | /(Module name)/(ID)/(Register name)/logstart()/(buffer size in bytes) Here: (buffer size in bytes) > 0 will turn on logging register changes (buffer size in bytes) = 0 will turn off logging register changes |
| Receive values with timestamp from the log | /(Module name)/(ID)/(Register name)/logget()/(number of records) |

This page is intentionally left blank

This chapter describes the setup of laser control over LAN or Wireless LAN.

Use the LAN connector on the front panel.

The WiFi module is optional. If applicable, the WiFi antenna is attached to the SMA connector labelled 'ANT'.

5.1. LAN Configuration

The server is initially configured for dynamic IP addressing. First, it is necessary to connect the laser to a local network running a DHCP server. After this, it is possible to change the configuration, once the communication has been established. Please note that a direct connection to PC via a LAN cable will not work, unless the PC is running DHCP server software.

1. Connect LAN connector to the Ethernet switch or router.
2. Connect PC to the same network.
3. Turn the laser ON.
4. Enter <http://Serial Number of Your Product> (e.g. <http://PGL426>) into the browser address bar.
5. The Main setup page is displayed (Figure 3).



The screenshot displays a web interface with several underlined navigation links in purple: Version, Set date and time, TCP config, REST, REMOTECONTROL.CSV, and RESTART. Below the REMOTECONTROL.CSV link, there is a file upload section containing a 'Browse...' button, the text 'REMOTECONTROL.CSV', a 'Submit' button, and a 'Reset' button.

Figure 3 Main setup page

6. Click *TCP config* to open the Network setup page.
7. Use the Network setup page to enable/disable LAN and WLAN interfaces and configure connection parameters.
8. Click *Submit* to save settings.

Done. Restarting... will appear on the screen.

After some time (usually tens of seconds), revisit the Main setup page by entering the URL discussed above.

If the wireless connection is used, add 'W' to the end of your URL, e.g. <http://PGL426W>.

5.2. Important Notes

- LAN connections are not password-protected.
- Any connection to the laser via LAN or WLAN in order to control it remotely must be done on a dedicated, limited access network (without access to the Internet). Connecting to a general office network will introduce hazards and risks to the device and personnel. This is especially important for wireless connections, as there is no easy way to quickly terminate communication (like unplugging a cable).
- The wireless module is radio equipment. Often there are import restrictions on radio equipment; therefore, the server is normally shipped without the wireless module installed.

"(-1) Log FIFO overrun"

"(-1) 3rd argument is missing (/s/s/???)"

"(-1) 2nd and 3rd arguments are missing (/s/???/???)"

"(1) End of registers list or enumerated values list"

"(2) No config file found"

"(3) wrong CFG file"

"(4) Application provided return buffer is too short"

"(5) No such device name"

"(6) No such register name"

"(7) Can't connect"

"(8) Timeout waiting for device answer"

"(9) Register is read only"

"(10) Register is not NV capable"

"(11) Violating top value limit"

"(12) Violating bottom value limit"

"(13) Wrong value, not included in allowed values list"

"(14) Register is not being logged"

"(15) Not enough memory"

"(16) Queue is empty"

"(17) Already connected, please disconnect first"

"(18) Not connected, please connect first"

This page is intentionally left blank

This chapter describes establishing, maintaining and troubleshooting communication with Ekspla lasers fitted with the Module.

It is assumed that installation and connecting procedures are accomplished.

7.1. Physical Connection

The Converter Module provides the following physical connectors on the front panel:

- USB DEVICE – the PC uses the USB communication devices class (CDC) driver included in Windows.
- RS232 – DCE, only Rx/D, Tx/D and ground are used.
- LAN connector – RJ45, 100Mb Ethernet.

At least one connection should be used, but pairing USB+LAN or RS232+LAN is recommended. In case LAN is not used, use USB, RS232, or both.

7.2. Testing and Troubleshooting Process Flow

1. Examination of the status of interfaces using the Windows Management Instrumentation (WMI) tool.
2. Probing of RS232 or virtual COM through USB connection.
3. Reading of data.

7.3. Communication Testing Utility

All three steps listed above are performed by the dedicated *CommTest.exe* utility. It uses the WMI service to discover COM ports, identifies the laser USB virtual COM port, probes the RS232 connection, scans the LAN for Ekspla lasers, etc.

The utility is located in the \Tools folder. It should be copied to your local hard disk drive and run from there.

After the utility is closed, a *log.txt* file is created in the same folder. The log file contains responses from WMI queries. The log file should be sent to Ekspla together with questions if support is needed.

7.4. Step-by-step Instructions

1. Connect the USB, LAN cables.
2. Run the *CommTest.exe*.

- Press the *Start inspection of hardware interfaces* button. If recognized the COM port will appear in *Working direct RS232 interfaces found* list.




Figure 4 *Working direct RS232 interfaces found* list

Please note that the entry *Ekspla RS232/CAN over USB device port* in the list is just a preliminary guess by the software. There is no reliable way to determine this using only WMI tools. Further steps describe the probing of ports.

If the serial driver is not loaded, a *Faulty interfaces found* message will appear.

When using Windows 10, the procedure should complete regardless of the serial driver loading status. For earlier Windows versions, loading of the serial driver depends on the success of the compatible driver search procedure. This procedure requires an internet connection and a fully updated Windows PC, with no updates/restarting pending or deferred. If the search procedure fails, install the *linux-cdc-acm.inf* file from

<https://www.kernel.org/doc/Documentation/usb/linux-cdc-acm.inf>

The .inf file instructs Windows to explicitly use *usbserial.sys* for the virtual COM port.

- Press *Next*. Working serial interfaces are listed. Choose the port or use the suggested one.
- Press *Probe Selected Port* to check communication. This sends */id()*, */list()* and */ip()* commands to the selected serial port. The responses are displayed.




Figure 5 *Successful connection response*

If there is no response to the `/ip()` command, try another COM port: select another one and press *Probe Selected Port* again. Try disconnecting and connecting the USB again.

For RS232 testing, check electrical wiring first. See connector pin signals on the laser side.




Figure 6 RS232 connector pinout on the laser side.

7.5. Testing LAN without an USB/RS232 Connection

1. Go to the *LAN* tab.
2. Press the *Get PC IP address(es)* button. This queries your PC's IP and fills the list. It is also possible to input IP addresses by hand.
3. Press the *Start LAN scan* button. It will attempt to connect to [`http://xxx.xxx.xxx.xxx:8080/id\(\)`](http://xxx.xxx.xxx.xxx:8080/id()) URLs throughout all the IP range.

The connection timeout is relatively short (100 ms) for the scan to be quick. This is usually enough for discovery in the LAN.

Get PC IP address(es)

192.168.0.51.fe80::f98c:870c:

Scan IP range

From:

192

168

0

100

To:

200

Start LAN scan

| Http adress | Device id() |
|----------------------------------|-------------------------------|
| http://192.168.0.194:8080/list() | SY+LD+PHD Date: 2015-12-28 |

Figure 7 LAN scan results

8.1. REST API over LAN or WLAN Register List Example

Table 6 REST API over LAN or WLAN register list example

| Name | ID | Reg ID | Menu | Type | User rights | Non-volatile | Min value | Max value | Short name | Print format | Name | Value |
|-------|------|--------|------|------|-------------|--------------|-----------|------------|------------|------------------------------------|----------------------------|---------------|
| LDD1A | \$12 | \$10 | \$0 | u8 | AUS | NV | 0 | 1 | State | [OFF,ON,FAULT] | Power | FAULT |
| LDD1A | \$12 | \$12 | \$2 | u16 | AUrS | NV | 0 | 2500 | Iset | %.3fA | Set Current | 0.850 |
| LDD1A | \$12 | \$47 | \$2 | u16 | AUrS | NV | 0 | 2500 | Istart | %.3fA | Set Start Current | 0.940 |
| LDD1A | \$12 | \$13 | \$3 | u16 | ArUrSr | | 0 | 1000 | Idisp | %.3fA | Display Current | 0.021 |
| LDD1A | \$12 | \$15 | \$2 | u16 | AUrS | NV | 0 | 1060 | Tstat | %.1fC | Thermostat | 100.0 |
| LDD1A | \$12 | \$19 | \$9 | u8 | ArUrSr | | 0 | 3 | Fault | [NONE, DRIVER, COOLING, INTERLOCK] | Fault source | COOLING |
| LDD1A | \$12 | \$1A | \$A | u16 | ArUrSr | | 0 | 65535 | Fcode | %xHEX | Fault code | 400 |
| LDD1A | \$12 | \$1B | \$B | u32 | AUrS | NV | 0 | 4294967295 | Time | %us | Work seconds | 22812090 |
| LDD1A | \$12 | \$28 | \$23 | u8 | ArUrSr | | 0 | 2 | IntLED | [OK,NOT OK,Defeated] | External interlock state | Defeated |
| LDD1A | \$12 | \$29 | \$24 | u8 | ArUrSr | | 0 | 2 | IntICD | [OK,NOT OK,Defeated] | Cooling interlock state | NOT OK |
| LDD1A | \$12 | \$2A | \$31 | u16 | ArUrSr | | 0 | 1023 | PhdOut | %u | Feedback photodiode output | 20 |
| LDD1A | \$12 | \$32 | \$41 | s16 | ArUrSr | | -2744 | 5000 | Tdisp | %.2fC | Display temperature | 47.38 |
| LDD1A | \$12 | \$34 | \$41 | u8 | ArUrSr | | 0 | 1 | Stable | [Not Stable,Stable] | Stable | Not Stable |
| LDD1A | \$12 | \$9 | \$F | u32 | ArUrSr | | 0 | 4294967295 | Version | %06u | Firmware | LDD1A; 131112 |
| LDD1A | \$12 | \$2 | \$0 | u16 | ArUrSr | | 0 | 65535 | ErrCode | %04x | Error Code | |

8.2. ASCII Serial Protocol Register List Example

Table 7 ASCII serial protocol register list example

| <i>Module name</i> | <i>Module ID</i> | <i>Type</i> | <i>User rights</i> | <i>Non-volatile</i> | <i>Min value</i> | <i>Max value</i> | <i>Print format</i> | <i>Register name</i> | <i>Captured value</i> | <i>Comments</i> |
|--------------------|------------------|-------------|--------------------|---------------------|------------------|------------------|-----------------------------|---|-----------------------|--|
| PHD1K000 | 48 | u16 | ArUrSr | | 0 | 65535 | %u | Data | 31956 | 1064nm output energy sensor, raw ADC data. Sensor is located before attenuator |
| PHD1K000 | 48 | float | ArUrSr | | 0 | 3,40 E+52 | %f | Mean | 100.997 | 1064nm energy sensor, calibrated and averaged for 1sec |
| SY3PL50M | 32 | u8 | AUS | NV | 0 | 1 | [OFF,ON, Failure] | State | ON | Laser state. OFF - laser stopped, ON - laser is firing. 'RUN' and 'STOP' buttons on the control pad change this register |
| SY3PL50M | 32 | u8 | AS | | 0 | 2 | [Continuous, Burst,Trigger] | Continuous / Burst mode / Trigger burst | Continuous | Burst mode control and trigger the burst |
| SY3PL50M | 32 | u32 | AUS | NV | 2 | 524287 | %u 1/OptClk | PRE-T delay | 2 | Allows preset the Pre-T OUT pulse delay |
| SY3PL50M | 32 | u16 | AUS | NV | 10 | 65535 | %.1fns | OUT3 delay | 14.Rgs | TRIG1 OUT pulse delay |
| SY3PL50M | 32 | u16 | AS | | 1 | 50000 | %u | Burst length, pulses | 1 | Allows setting the number of pulses in the burst to be emitted when operating laser in burst mode |
| SY3PL50M | 32 | u8 | AUrS | NV | 0 | 1 | [Internal, External] | Synchronization Mode | Internal | When External mode is set, laser is triggered from external signal, supplied to SYNC IN |
| SY3PL50M | 32 | u8 | AUS | NV | 0 | 2 | [OFF, Adjustment, Maximum] | Energy level | OFF | Operating mode. See 'Parameters Description' in 5 Chapter of Manual for description |

| <i>Module name</i> | <i>Module ID</i> | <i>Type</i> | <i>User rights</i> | <i>Non-volatile</i> | <i>Min value</i> | <i>Max value</i> | <i>Print format</i> | <i>Register name</i> | <i>Captured value</i> | <i>Comments</i> |
|--------------------|------------------|-------------|--------------------|---------------------|------------------|------------------|-------------------------|---------------------------|-----------------------|--|
| SY3PL50M | 32 | u16 | AS | NV | 1 | 5000 | %u | Frequency divider | 1 | Laser pulse repetition rate is divided by a factor set in this menu. Pumping frequency remains unchanged |
| SY3PL50M | 32 | u16 | AUS | NV | 50 | 5000 | %uus | Pump delay, adj. level | 348 | Delay of regen pump in adjustment mode. Allows tweak pulse energy in adjustment mode |
| SY3PL50M | 32 | u32 | ArUrSr | | 0 | 4294967295 | %uHz | Optical Clock | 87551104 | Oscillator frequency, |
| SY3PL50M | 32 | u16 | ArUrSr | | 0 | 1000 | %.1fHz | External Syncln frequency | 1000.1 | Repetition rate, measured |
| SM5 | 61 | s32 | ArUrSr | | -2,00 E+09 | 2147483647 | %d | Current position | 261 | Attenuator stepper motor actual position |
| SM5 | 61 | s32 | AUS | NV | -2,00 E+09 | 2147483647 | %d | Target position | 261 | Attenuator stepper motor target position |
| SM5 | 61 | u8 | AUS | NV | 0 | 1 | [POSITIONING, VELOCITY] | Mode | POSITIONING | Attenuator stepper motor driver mode |
| CPU8000 | 17 | u16 | ArUrSr | | 0 | 1300 | %.1fA | Display Current | 0.4 | Laser diagnostic register |
| HV40W | 40 | u16 | ArUrSr | | 0 | 4000 | %uV | HV voltage | 11 | Laser diagnostic register |
| LDCO48BP | 28 | s16 | ArUrSr | | -200 | 4600 | %.2fC | Display temperature | 28.64 | Laser diagnostic register |
| LDM6A | 16 | u16 | ArUrSr | | 0 | 1000 | %.2fA | Display Current | 1.97 | Laser diagnostic register |
| LDM6A | 16 | s16 | ArUrSr | | -2744 | 5000 | %.2fC | Display temperature | 30.98 | Laser diagnostic register |

8.3. List of Commands for REST over API Protocol

Table 8 List of commands for REST over API protocol

| <i>Communication test</i> | |
|---------------------------|------------------|
| Command | / [slash symbol] |

| | |
|----------------------------|---|
| Sample response | Remote control REST app (Nov 10 2015). |
| HTML format | <html><TABLE> <tr><th>Remote control REST app (Nov 10 2015). </th></tr> </TABLE></html> |
| Device ID query | |
| Command | /id() Sends the first line from REMOTECONTROL.CSV file. |
| Sample response | Device ID SY320100 Date: 2015.10.29 Here, 'SY320100' is the device type, 'Date' is the configuration file creation date. |
| HTML format | <html><TABLE> <tr><th>Device ID</th></tr> <tr><td id=DN>SY320100 Date: 2015.10.29</td></tr> </TABLE></html> |
| Register list query | |
| Command | /list() Sends list of registers from REMOTECONTROL.CSV file, together with their respective URLs. |
| Sample response | IO:15 Error Code RA0,pin22,AN0 RA0 direction RA1,pin23,AN1 !SY3PL50M:32 Error Code State Synchronization Mode External Syncln frequency Photodiode bias Optical Clock |
| HTML format | <TABLE> <tr><th>IO:15</th></tr> <tr><td>Error Code</td></tr> <tr><td>RA0,pin22,AN0</td></tr> <tr><td>RA0 direction</td></tr> <tr></tr> <tr><th>SY3PL50M:32</th></tr> <tr><td>Error Code</td></tr> <tr><td>State</td></tr> <tr><td>Synchronization Mode</td></tr> <tr><td>External Syncln frequency</td></tr> <tr><td>Photodiode bias</td></tr> <tr><td>Optical Clock</td></tr> <tr></tr> </TABLE> |

| <i>Read command (numerical type variable)</i> | |
|--|--|
| Command | /(Module name)/(ID)/(Register name) E.g. /SY3PL50M/32/Oscillator clock |
| Sample response | Get register Device SY320100:32 Register Optical Clock Min. value 0 Max. value 4.29497e+09 RW No NV No Format %uHz Error (0) Success, no error Value 0Hz |
| HTML format | <pre> <html><body><TABLE> <tr><th>Get register</th></tr> <tr><th></th><th></th></tr> <tr><td>Device</td><td id=D1>SY320100:32</td></tr> <tr><td>Register</td><td id=R1>Optical Clock</td></tr> <tr><td>Min. value</td><td id=minV>0</td></tr> <tr><td>Max. value</td><td id=maxV>4.29497e+09</td></tr> <tr><td>RW</td><td id=RW>No</td></tr> <tr><td>NV</td><td id=Nv>No</td></tr> <tr><td>Format</td><td id=FMT>%uHz</td></tr> <tr><td>Error</td><td id=E1>(0) Success, no error</td></tr> <tr><td>Value</td><td id=V1>0Hz</td></tr></TABLE></body></html> </pre> |
| <i>Read command (set type variable)</i> | |
| Command | /(Module name)/(ID)/(Register name) E.g. /SY3PL50M/32/State |

| | |
|-----------------------------------|--|
| Sample response | Get register Device SY320100:32 Register Command Min. value 0 Max. value 4 RW Yes NV Yes Format [POWEROFF,SLEEP,STOP,PAUSE,RUN,FAULT] Error (0) Success, no error Value FAULT |
| HTML format | <pre> <html><body><TABLE><tr><th>Get register</th></tr> <tr><th></th><th></th></tr> <tr><td>Device</td><td id=D1>SY320100:32</td></tr> <tr><td>Register</td><td id=R1>Command</td></tr> <tr><td>Min. value</td><td id=minV>0</td></tr> <tr><td>Max. value</td><td id=maxV>4</td></tr> <tr><td>RW</td><td id=RW>Yes</td></tr> <tr><td>NV</td><td id=NV>Yes</td></tr> <tr><td>Format</td><td id=FMT>[POWEROFF,SLEEP,STOP,PAUSE,RUN,FAULT]</td></tr> <tr><td>Error</td><td id=E1>(0) Success, no error</td></tr> <tr><td>Value</td><td id=V1>FAULT</td></tr></TABLE></body></html> </pre> |
| Write command | |
| Command | (Module name)/(ID)/(Register name)/(Value) E.g. /SY320100/32/Repetition rate/500 This example command will set the oscillator clock to 500 kHz. |
| Sample response, accepted command | Set register to Device SY320100:32 Register Repetition rate Value 500 Error (0) Success, no error |
| HTML format, accepted command | <pre> <body><TABLE><tr><th>Set register to</th></tr> <tr><th></th><th></th></tr> <tr><td>Device</td><td id=D1>SY320100:32</td> <tr><td>Register</td><td id=R1>Repetition rate</td> <tr><td>Value</td><td id=V1>500</td> <tr><td>Error</td><td id=E1>(0) Success, no error</td></tr></TABLE></body> </pre> |

| | |
|--|---|
| Sample response, rejected command | Set register to Device SY320100:32 Register Repetition rate Value 1005 Error (11) Violating top value limit |
| HTML format, rejected command | <pre> <body><TABLE> <tr><th>Set register to</th></tr> <tr><th></th><th></th></tr><tr><td>Device</td><td id=D1>SY320100:32</td> <tr><td>Register</td><td id=R1>Repetition rate</td><tr><td>Value</td><td id=V1>1005</td> <tr><td>Error</td><td id=E1>(11) Violating top value limit</td> </tr></TABLE></body> </pre> |
| Write to NVRAM command format | |
| Command | /(Module name)/(ID)/(Register name)/(Value)/NV E.g. /SY320100/32/Repetition rate/500/NV |
| Sample response | Set NV register to Device SY320100:32 Register Repetition rate Value 500 Error (0) Success, no error |
| HTML format | <pre> <body><TABLE><tr><th>Set NV register to</th></tr> <tr><th></th><th></th></tr> <tr><td>Device</td><td id=D1>SY320100:32</td> <tr><td>Register</td><td id=R1>Repetition rate</td> <tr><td>Value</td><td id=V1>500</td> <tr><td>Error</td><td id=E1>(0) Success, no error</td> </tr></TABLE></body> </pre> |
| Start/stop logging register updates | |
| Command | /(Module name)/(ID)/(Register name)/logstart()/(buffer size in bytes) E.g. /SY320100/32/Photodiode%20bias/logstart()/100 (buffer size in bytes) > 0 will turn on logging of register changes (buffer size in bytes) = 0 will turn off logging of register changes The sample command will turn on logging of register 'Photodiode bias'. |
| Sample response, accepted command | Start register log Device SY320100:32 Register Photodiode bias Log memory size 100 Error (0) Success, no error |

| | |
|---|---|
| HTML format, accepted command | <pre> <body><TABLE><tr><th>Start register log</th></tr> <tr><th></th><th></th></tr> <tr><td>Device</td><td id=D1>SY320100:32</td> <tr><td>Register</td><td id=R1>Photodiode bias</td> <tr><td>Log memory size</td><td id=V1>100</td> <tr><td>Error</td><td id=E1>(0) Success, no error</td> </tr></TABLE></body> </pre> |
| Receive values with timestamp from the log | |
| Command | <pre> /(Module name)/(ID)/(Register name)/logget()/(number of records) E.g. /SY3PL50M/32/Photodiode bias/logget()/5 </pre> |
| Sample response, <u>buffer overrun error</u> | <pre> Get register log Device SY320100:32 Register Photodiode bias Error Log overflow detected Data Timestamp [ms] 5.93V 530169612 5.93V 530169633 5.93V 530169653 5.93V 530169674 5.93V 530169695 </pre> |
| HTML format, <u>buffer overrun error</u> | <pre> <body><TABLE><tr> <th>Get register log</th> </tr> <tr> <th></th> <th></th> </tr> <tr> <td>Device</td> <td id=D1>SY320100:32</td> </tr> <tr> <td>Register</td> <td id=R1>Photodiode bias</td> </tr> <tr> <td>Error</td> <td id=E1>Log overflow detected</td> </tr> <tr> </tr></TABLE><TABLE id=L1><tr> <th>Data</th> <th>Timestamp [ms]</th> </tr> <tr> <td>5.93V</td> <td>530169612</td> </tr> <tr> <td>5.93V</td> <td>530169633</td> </tr> <tr> <td>5.93V</td> <td>530169653</td> </tr> <tr> <td>5.93V</td> <td>530169674</td> </tr> <tr> <td>5.93V</td> <td>530169695</td> </tr> <tr> <td></td> <td></td> </tr><tr> <td></td> <td></td> </tr> <tr> <td></td> <td></td> </tr><tr> <td></td> <td></td> </tr> <tr> <td></td> <td></td> </tr></TABLE></body> </pre> |
| Command | <pre> /(Module name)/(ID)/(register name)/logget()/(number of records) E.g. /CAMERA/57/Mean/logget()/10 </pre> |

| | |
|-----------------------------------|--|
| Sample response, <u>no errors</u> | Get register log Device CAMERA:57 Register Mean Error (0) Success, no error Data Timestamp [ms] 384.000000 535959608 384.000000 535960119 384.000000 535960630 384.000000 535961142 384.000000 535961653 |
| HTML format, <u>no errors</u> | <pre> <TABLE><tr> <th>Get register log</th> </tr> <tr> <th></th> <th></th> </tr><tr> <td>Device</td> <td id=D1>CAMERA:57</td> </tr> <tr> <td>Register</td> <td id=R1>Mean</td> </tr> <tr> <td>Error</td> <td id=E1>(0) Success, no error</td> </tr> <tr> </tr></TABLE><TABLE id=L1><tr> <th>Data</th> <th>Timestamp [ms]</th> </tr> <tr> <td>384.000000</td> <td>535959608</td> </tr> <tr> <td>384.000000</td> <td>535960119</td> </tr> <tr> <td>384.000000</td> <td>535960630</td> </tr> <tr> <td>384.000000</td> <td>535961142</td> </tr> <tr> <td>384.000000</td> <td>535961653</td> </tr> <tr> <td></td> <td></td> </tr><tr> <td></td> <td></td> </tr> <tr> <td></td> <td></td> </tr><tr> <td></td> <td></td> </tr> <tr> <td></td> <td></td> </tr></TABLE> </pre> |

8.4. List of Commands for ASCII Serial Protocol

Table 9 List of commands for ASCII serial protocol

| Communication test | |
|--------------------|--|
| Command | [CR] |
| Sample response | Remote control over RS232 (Jun 18 2015) Here, 'Jun 18 2015' is the version date of the RS232 interpreter. |
| Device ID query | |
| Command | /id() [CR] |

| | | |
|-----------------------------------|---|--|
| Sample response | Device: DNL207 Date: 17/09/2015 Here, 'DNL207' is the device serial number, 'Date' is the configuration file creation date. | |
| Register list query | | |
| Command | /list() [CR] | |
| Sample response | LDM6A:16 Error Code Display Current Display temperature UCP:4 Error Code Firmware UCPv3 HV40W:40 Error Code HV voltage TK6:44 Error Code Display temperature SM5:61 Error Code Current position Target position Mode CPU8000:17 Error Code Display Current LDCO48BP:28 Error Code Display temperature SY3PL50M:32 Error Code State Continuous / Burst mode / Trigger burst PRE-T clock divider PRE-T delay OUT3 delay Burst length, pulses Synchronization Mode Energy level Frequency divider Output redirection Pump delay, adj. level Optical Clock External SyncIn frequency PHD1K000:48 Error Code Data Mean Batch size PS5050:22 Error Code Channel 1 state | |
| Read command | | |
| Command | /(Module name)/(ID)/(Register name) [CR] E.g. /SY3PL50M/32/State | |
| Sample response | ON#13#10#03 | |
| Write command | | |
| Command | (Module name)/(ID)/(Register name)/(Value) [CR] E.g. /SY3PL50M/32/State/ON This example command will turn on the laser. | |
| Sample response, accepted command | #13#10#03 | |

| | |
|--|---|
| Sample response, rejected command | ""Error: (13) Wrong value, not included in allowed values list#13#10#03 |
| Write to NVRAM command format | |
| Command | /(Module name)/(ID)/(Register name)/(Value)/NV [CR] E.g. /SY3PL50M/32/State/ON/NV This example command will turn the laser on and will make the ON value non-volatile, i.e. at the next power-on, the laser will try to run automatically. |
| Sample response, accepted command | #13#10#03 |
| Sample response, rejected command | ""Error: (13) Wrong value, not included in allowed values list#13#10#03 |
| Start/stop logging register updates | |
| Command | /(Module name)/(ID)/(Register name)/logstart()/(buffer size in bytes) [CR] E.g. /CPU8000/17/Display Current/logstart()/2000 (buffer size in bytes) > 0 will turn on logging of register changes (buffer size in bytes) = 0 will turn off logging of register changes The sample command will turn on logging of register 'Display Current'. |
| Sample response, accepted command | #13#10#03 |
| Receive values with timestamp from the log | |
| Command | /(Module name)/(ID)/(Register name)/logget()/(number of records) [CR] E.g. /CPU8000/17/Display Current/logget()/100 |
| Sample response, no errors, messages present in the buffer | (value1) (time stamp1)#13#10 (value2) (time stamp2)#13#10 (value3) (time stamp3)#13#10 (value4) (time stamp4)#13#10 #03 0.020A 634698780#13#10 0.020A 634698799#13#10 0.020A 634698821#13#10 0.020A 634698840#13#10 |

| | |
|-----------------------------------|---|
| Sample response, empty buffer | #13#10#03 |
| Sample response, buffer overrun | 0.020A 634698780#13#10 0.020A 634698799#13#10 0.020A 634698821#13#10 0.020A 634698840#13#10 0.020A 634698862#13#10 0.020A 634698881#13#10 ""Error: (-1) Log FIFO overrun#13#10#03 |
| Sample response, rejected command | ""Error: (14) Register is not being logged#13#10#03 |
| Device timer value query | |
| Command | /timestamp [CR] |
| Sample response | Timer value#13#10#03 Here, 'Timer value' is an actual value from the timer, e.g. 88479777. |